

TABLE 9.1 Examples of Cell Immobilization by Entrapment Using Different Support Materials

Cells	Support matrix	Conversion
<i>S. cerevisiae</i>	κ -Carrageenan or polyacrylamide	Glucose to ethanol
<i>E. aerogenes</i>	κ -Carrageenan	Glucose to 2,3-butanediol
<i>E. coli</i>	κ -Carrageenan	Fumaric acid to aspartic acid
<i>Trichoderma reesei</i>	κ -Carrageenan	Cellulose production
<i>Z. mobilis</i>	Ca-alginate	Glucose to ethanol
<i>Acetobacter</i> sp.	Ca-alginate	Glucose to gluconic acid
<i>Morinda citrifolia</i>	Ca-alginate	Anthraquinone formation
<i>Candida tropicalis</i>	Ca-alginate	Phenol degradation
<i>Nocardia rhodocrous</i>	Polyurethane	Conversion of testosterone
<i>E. coli</i>	Polyurethane	Penicillin G to G-APA
<i>Catharantus roseus</i>	Polyurethane	Isocitrate dehydrogenase activity
<i>Rhodotorula minuta</i>	Polyurethane	Menthyl succinate to menthol

In mixed-culture microbial films, the presence of some polymer-producing organisms facilitates biofilm formation and enhances the stability of the biofilms. Microenvironmental conditions inside a thick biofilm vary with position and affect the physiology of the cells.

In a stagnant biological film, nutrients diffuse into the biofilm and products diffuse out into liquid nutrient medium. Nutrient and product profiles within the biofilm are important factors affecting cellular physiology and metabolism. A schematic of a biofilm is depicted in Fig. 9.11. Biofilm cultures have almost the same advantages as those of the immobilized cell systems over suspension cultures, as listed in the previous section.

The thickness of a biofilm is an important factor affecting the performance of the biotic phase. Thin biofilms will have low rates of conversion due to low biomass concentration, and thick biofilms may experience diffusionally limited growth, which may or may not be beneficial depending on the cellular system and objectives. Nutrient-depleted regions may also develop within the biofilm for thick biofilms. In many cases, an optimal biofilm thickness resulting in the maximum rate of bioconversion exists and can be determined. In some cases, growth under diffusion limitations may result in higher yields of products as a result of changes in cell physiology and cell-cell interactions. In this case,

TABLE 9.2 Examples of Cell Immobilization by Surface Attachment

Cells	Support surface	Conversion
<i>Lactobacillus</i> sp.	Gelatin (adsorption)	Glucose to lactic acid
<i>Clostridium acetobutylicum</i>	Ion-exchange resins	Glucose to acetone, butanol
<i>Streptomyces</i>	Sephadex (adsorption)	Streptomycin
Animal cells	DEAE-sephadex/cytodex (adsorption)	Hormones
<i>E. coli</i>	Ti(IV) oxide (covalent binding)	
<i>B. subtilis</i>	Agarose-carbodiimide (covalent binding)	
<i>Solanum aviculare</i>	Polyphenylene oxide-glutaraldehyde (covalent binding)	Steroid glycoalkaloids formation